

UPRIGHT REFRIGERATORS AND FREEZERS



Models:

ACR82L ACR1818 ACRR432L ACFF436L ACR261RH AFG26MLLH SCR801G SCR1802G SCRR432 SCFF436 SCRR261G SCFF262GLH

User Manual

BEFORE USE, PLEASE READ AND FOLLOW ALL SAFETY RULES AND OPERATING INSTRUCTIONS

Write Model and Serial Numbers here:
Model
Serial No

Felix Storch, Inc. An ISO 9001:2015 registered company 770 Garrison Ave Bronx, New York 10474 www.accucold.com

TABLE OF CONTENTS

IMPORTANT SAFEGUARDS	3
INSTALLATION INSTRUCTIONS5	
Before Using Your Appliance5	5
Locating your Appliance5	5
Installing the Leveling Legs6	
Installing the Shelves7	7
Installing the Pressure Valve7	
Connecting the Power Supply7	7
Door Reversal8	3
OPERATING YOUR APPLIANCE 1	1
CARE AND MAINTENANCE 1	4
TROUBLESHOOTING CHART 1	5

<u>CALIFORNIA CARB/SNAP DISCLOSURE</u> This product uses eco-friendly hydrocarbon refrigerant and fully complies with California CARB regulations.

However, we are required by California Law to provide the following disclosure statement in every product sold in California.

"This equipment is prohibited from use in California with any refrigerants on the 'List of Prohibited Substances' for that specific end-use, in accordance with California Code of Regulations, title 17, section 95374. This disclosure statement has been reviewed and approved by Felix Storch, Inc. and Felix Storch, Inc. attests, under penalty of perjury, that these statements are true and accurate."

"This product does <u>not</u> use any refrigerants on the 'List of Prohibited Substances'"

IMPORTANT SAFEGUARDS

Your safety and the safety of others are very important.

We have attached many important safety messages in this manual and on your appliance. Always read and obey all safety messages.



This is the Safety Alert Symbol. This symbol alerts you to potential hazards that can kill or injure you and others. All safety messages will follow the Safety Alert Symbol and either the words "DANGER" "WARNING" or "CAUTION".

WARNING

DANGER means that failure to heed this safety

statement may result in severe personal injury or death.

DANGER

WARNING means that failure to heed this safety statement may result in extensive product damage, serious personal injury, or death.

CAUTION

CAUTION means that failure to heed this safety statement may result in minor or moderate personal injury, or property or equipment damage.

Safety messages will alert you to what the potential hazard is, tell you how to reduce the chance of injury, and let you know what can happen if the user does not follow the instructions.



Before using your refrigerator/freezer, it must be positioned properly and installed as described in this manual, so read the manual carefully. To reduce the risk of fire, electrical shock, or injury when using the refrigerator/freezer, follow basic precautions, including the following:

DANGER

- These appliances are CFC and HFC free and contain a small quantity of liquid propane (R290) which is environmentally friendly, but flammable. It does not damage the ozone layer or increase the greenhouse effect. Care must be taken during transportation and while setting up the appliance so that no cooling system parts are damaged. Leaking refrigerant is combustible and may damage the eyes.
- In case of any damage:
 - o Avoid open flames and anything that creates a spark
 - o Disconnect from the power source
 - o Ventilate the room in which the appliance is in for several minutes and contact the Service Department for advice
- The more refrigerant in an appliance, the larger the room it is being installed in should be. In case of a leakage, if the appliance is in a small room, there is a greater risk of combustible gases building up. The amount of refrigerant in the appliance can be found on the rating label on the appliance. It is hazardous for anyone other than the Authorized Service Person to conduct servicing or repairs to the appliance.
- Take serious care when handling, moving, and using the appliance to avoid damaging the refrigerant tubing and increasing the risk of a leak.
- Parts replacement and servicing shall be done by factory authorized service personnel to minimize the risk of ignition due to incorrect parts or improper service.

- Plug into a grounded 3-prong outlet. Do not remove the grounding prong, do not use an adapter, and do not use an extension cord.
- Replace all panels before using the unit.
- lt is recommended that a separate circuit, serving only your refrigerator/freezer, be provided. Use an electrical outlet that cannot be turned off by a switch or pull chain.
- Never clean refrigerator/freezer parts with flammable fluids, these fumes can create a fire hazard or explosion.
- Do not store or use gasoline or other flammable vapors and liquids in the vicinity of this or any other appliance. The fumes can create a fire hazard or explosion.
- Before starting cleaning and maintenance operations, make sure the unit's power line is disconnected. Do not connect or disconnect the electric plug when your hands are wet.
- Make sure to disconnect the unit's power supply before starting cleaning and maintenance operations, failure to do so can result in electrical shock or death.
- Do not attempt to repair or replace any part of your refrigerator/freezer unless it is specifically recommended in this material. A qualified technician should complete all other servicing.

WARNING

This unit is not designed to operate in enclosed spaces. When placing your unit, make sure you allow adequate airflow on all sides.

- Use two or more people to move and install the refrigerator/freezer. Failure to do so can result in back or other injuries.
- To ensure proper ventilation for your refrigerator/freezer, the front of the unit must be completely unobstructed. Choose a well-ventilated area with temperatures above 59°F (15°C) and below 86°F (30°C). This unit must be installed in an area protected from the elements, such as wind, rain, water spray, or sunlight.
- The refrigerator/freezer should not be installed next to ovens, grills, or other sources of high heat.
- The refrigerator/freezer must be installed with all electrical connections in accordance with state and local codes. This unit uses the standard electrical supply (115 V AC, 60hz) and must be properly grounded in accordance with national and local codes and ordinances.
- Do not kink or pinch the power supply cord of the refrigerator/freezer.
- The fuse (or circuit breaker) size should be 15 amperes.
- It is important for the refrigerator/freezer to be leveled in order for it to function properly. You may need to make several adjustments to level it.
- Never allow children to use, play with or crawl inside the refrigerator/freezer
- Do not use solvent-based cleaning agents or abrasives on the interior, these cleaners may damage or discolor the interior.
- Keep your fingers out of the "pinch point" areas. Clearances between the door and cabinets are necessarily small. Be careful closing the door(s) when children are in the area.
- Do not use this appliance for anything other than its intended purpose.



Risk of child entrapment

Child entrapment and suffocation are not problems of the past. Junked or abandoned freezers are still dangerous, even if they will "just sit in the garage a few days."

- Before you throw away your old freezer: Take off the doors. Leave the shelves in place so that children may not easily climb inside.
- Never allow children to use, play with, or crawl inside the freezer.

- SAVE THESE INSTRUCTIONS -

INSTALLATION INSTRUCTIONS

Before Using Your Appliance

- Remove the exterior and interior packing.
- Check to be sure you have all the following parts:
 - Set of pop-out keys
 - 1 instruction manual
- Before connecting the appliance to the power source, let it stand upright for approximately 2 hours or let the compressor complete three on/off cycles to assure that the temperature inside is appropriate for the products before loading them into the appliance. This will reduce the possibility of a malfunction in the cooling system from handling during transportation.
- Clean the interior surface with lukewarm water using a soft cloth.

Locating your Appliance

- Move the appliance as close as possible to its final location, before removing the skid.
- Use a screwdriver to remove the screws that fix the skid to the cooler or metal supports, if applicable.
- Place the cabinet on a surface solid enough to prevent vibration and strong enough to support the combined weight of the cabinet and the product load. Estimate 35 pounds per cubic foot of storage space.
- Cabinets should never be placed in front of windows or glass doors where sunlight will be in direct contact with the units. If a cabinet is located near a window or glass door, provide an adequate shade to block the sun's rays.
- Ambient temperature must be between 59°F (15°C) and 86°F (30°C). Low ambient temperatures below 50° F are as harmful as high ambient temperatures. Extreme temperatures will affect the performance of the cabinet.
- **NOTE:** Failure to properly level the cabinet could hinder proper drainage and door operation. Leveling legs are provided with the cabinets. Level the cabinet by adjusting the leveling legs as needed.
- Leave a space of at least 8in (203mm) between the back, sides, and top of the cabinet to provide adequate air circulation of the unit compartment.

Installing the Leveling Legs

Heavy-duty leveling legs are included so the unit may be leveled accurately and to allow room to clean beneath the cabinet.

- CAUTION: DO NOT attempt to install the level legs by yourself. You are far less likely to incur injury to yourself or the appliance if you install the legs with another person.
- Installing caster wheels: caster wheels are included with some models. Only use genuine caster wheels or the base of the unit might suffer damage. When installing the caster wheels, it is important to tighten them against the base rail (see image) The top edge of the caster wheels must make firm contact with the base rail. The top edge of the caster wheels must be in firm contact with the base rail.



NOTE: The caster wheels with handbrakes must be placed in FRONT of the unit.

- 2. Heavy-duty leveling legs are provided so the unit can be leveled accurately and to allow room to clean beneath the cabinet.
- 3. Tilt the cabinet back just far enough to fit the front legs and place a block of wood or some similar support under it until both front legs are in. Then tilt the cabinet forward enough to install the rear legs.
- 4. After the legs are installed, level the cabinet by using the adjustable part of the legs as needed. This is very important to ensure proper draining and door operation.



Installing the Shelves

Shelves and shelf clips (4 clips per shelf) are packed inside the cabinet before shipment. After choosing the appropriate spacing inside the cabinet, install shelf clips by inserting the top of the clip inside the pilaster and then press the lower section of the clip while pushing inwards until it snaps into place. Repeat until all clips are in the desired location. Place the shelf on clips, checking that all four corners of the shelf are leveled and supported by all four shelf clips.



Installing the Pressure Valve

- 1. Locate the pressure relief valve packed inside of your appliance.
- 2. Locate the hole on the backside. (See image, right.)
- 3. Align the pressure relief valve with the hole, using the marking as a guide.
- 4. Push the pressure relief valve until you hear a click, this means the valve was correctly installed.

Connecting the Power Supply

Each unit requires a dedicated circuit to prevent the malfunction of the compressor and the electrical components. Each unit should be connected to its own circuit in accordance with National Electric Codes and local ordinances.



NOTE: DO NOT connect the appliance to a wall outlet with other appliances also plugged in. DO NOT use extension cords to connect the appliance to the power outlet.

Plug-in electrical power cord: Do not under any circumstance cut or remove the grounding pin from the power cord. Removal of grounding pin invalidates product warranty.

Door Reversal

APPLY ONLY FOR SINGLE DOOR MODELS (Except ACR82L/SCR801G)

TOOLS

If you would like to reverse the door swing, follow the steps below.

Thirester

5

3/8" NUT DRIVER BIT

PHILLIPS # 2 SCREWDRIVER

INSTRUCTIONS





Step 2: Remove the acrylic header sign and place it on the top of the unit.



Step 3: Remove the four screws from the header sign frame.



Step 4: Remove the four screws from the front grill, then remove the front grill from the cooler.



Step 5: Remove the screws of the upper and lower hinge. Hold the door when you are removing the screws.



Step 6: With the hinges installed on the door, rotate the door 180 degrees.



Step 7: Install the door on the cooler, header sign, and front grill according to steps 1 to 5.

OPERATING YOUR APPLIANCE

Control Panel

Your unit is manufactured with a Dixell electronic thermostat to control various functions of your unit, such as temperature, defrost cycles, alarms ,etc.



TEMPERATURE SELECTION

TYPE	TEMPERATUR E RANGE	DEFAUL T SET POINT	CONTROLLE R MODEL
REFRIGERATO R	+2°C to +8°C (+36°F to +46°F)	+5°C (+41°F)	DIXELL XR40CX
FREEZER	-20.5°C to -15°C (-5°F to +5°F)	-20°C (- 4°F)	DIXELL XR60CX

FRONT PANEL COMMANDS

KEY/BU	JTTON	FUNCTION
SET		To display the target set point. In programming mode, it selects a parameter or confirms an operation
444	DEF	To start a manual defrost
Δ	UP	To see the maximum stored temperature. In programming mode, it browses the parameter codes or increases the displayed value
\triangleleft	DOWN	To see the minimum stored temperature. In programming mode, it browses the parameter codes or decreases the displayed value.
\bigcirc		To turn the unit off
÷.		Not enabled

DISPLAY INDICATORS

LED MODE			FUNCTION
	*	ON	Compressor enabled
	*	Flashing	Anti-short cycle delay enabled
	****	ON	Defrost enabled
	*** ****	Flashing	* Drip time in progress
5		ON	* Fans enabled
5		Flashing	* Fans delay after defrosting in progress
«)		ON	* An alarm is occurring

* Freezers only

KEY COMBINATIONS:

$\triangle + \bigtriangledown$	To lock & unlock the keyboard	
SET +	To enter programming mode	
SET + A	To return to the room temperature display	

HOW TO SEE THE SET POINT:

- 1. The set point parameter is the temperature at which the compressor cuts out.
- 2. Push and immediately release the SET key: the display will show the set point value.
- 3. Push and immediately release the SET key or wait for 5 seconds to display the probe value again.

HOW TO CHANGE THE SET POINT:

- 1. Push the SET key for more than 5 seconds to change the SET point value.
- The value of the set point will be displayed and the "C" or "F" LED starts blinking.
- 3. To change the Set value,' push the \land or \checkmark arrows within 10 seconds.
- 4. To memorize the new set point value, push the SET key again or wait for 10 seconds.

HOW TO SEE THE MAXIMUM TEMPERATURE RECORDING

- 1. Press and release the \blacktriangle (UP) key.
- 2. The "Hi" message will be displayed followed by the maximum temperature recorded.
- 3. By pressing the o key again or by waiting 5s the normal display will be restored.

HOW TO SEE THE MINIMUM TEMPERATURE RECORDING

- 1. Press and release the ▼ (DOWN) key.
- 2. The "Hi" message will be displayed followed by the maximum temperature recorded.
- 3. By pressing the o key again or by waiting 5s the normal display will be restored.

HOW TO RESET THE MAX AND MIN TEMPERATURE RECORDED

- 1. Hold and press the SET key for more than 3s, while the max. or min temperature is displayed. (rSt message will be displayed)
- 2. To confirm the operation the "rSt" message starts blinking and the normal temperature will be displayed

			Free	ezer	Refrig	erator
Group	Parameter	Description	°C	°F	°C	°F
Regulation	Set	Set point	-20	-4	5	41
Regulation	Ну	Differential	1	2	0.5	1
Regulation	LS	Minimum set point	-20.5	-5	2	36
Regulation	US	Maximum set point	-15	5	8	46
Defrost	dtE	Defrost termination temperature first defrost	13	55	-	-
Fan	FSt	Fan stop temperature	37	99	37	99
Alarm	ALU	Maximum temperature alarm	0	32	15	59
Alarm	ALL	Minimum temperature alarm	-25	-13	0	32
Alarm	AFH	Differencial for temperature alarm recovery	1	2	0.5	1
Display	CF*	Temperature measurement unit				

HOW TO CHANGE A PARAMETER VALUE

The hidden menu includes all the parameters of the instrument.

HOW TO ENTER THE HIDDEN MENU

- 1. Enter the Programming mode by pressing the **SET+DOWN** keys for 3 seconds (the "°**C**" or "°**F**" LED starts blinking).
- 2. Release the keys, then again push the **SET+DOWN** keys for more than 7 seconds. The Pr2 label will be displayed immediately following the HY parameter.

NOW YOU ARE IN THE HIDDEN MENU.

- 3. Select the required parameter.
- 4. Press the "**SET**" key to display its value
- 5. Use UP or DOWN to change its value.
- 6. Press "SET" to store the new value and move to the following parameter.

To exit: Press SET+UP or wait for 15 seconds without pressing a key.

NOTE 1: If no parameter is present in Pr1, after 3 seconds the "noP" message is displayed, keep the keys pushed until the Pr2 message is displayed.

NOTE 2: the set value is stored even when the procedure is halted by waiting for the time-out to expire.

* **WARNING:** When the measurement unit is changed the SET point and the values of the following parameters Hy, LS, US, OT, ALU, and ALL must be checked and modified if necessary.

DIXELL ALARM INDICATORS

MESSAGE	CAUSE	COMMENTS	
НА	Maximum	This alarm is DEFAULT VALU	
	temperature alarm	adjustable. To change the values, it is	Refrigerator 15 °C (59 °F)
		necessary to enter the hidden menu and change the parameters.	Freezer 0°C (32 °F)
LA			DEFAULT VALUES
	temperature alarm adjustable. To change the values, it is necessary to enter the hidden menu and change the parameters.	Refrigerator 0°C (32 °F)	
		the hidden menu and change the	Freezer -25°C (-13 °F)
DA	Door alarm	If the door stays open for more than 3 minutes, an audible and visual alarm will go off on the thermostat	

CARE AND MAINTENANCE

Your unit requires little maintenance. However, periodic cleaning can help prolong product life. Remember to disconnect the power source before cleaning. Remove all products and place them in another cooler.

- 1. Clean the inside and the exterior of your cooler, including the shelves, with a mild soap solution. Do not use strong detergents. Wipe dry the inside of the cooler before loading products again.
- 2. Plug in your unit and wait 2 hours for the temperature to reach the appropriate range to store your products in.

TROUBLESHOOTING CHART

You can solve many common freezer problems easily, saving you the cost of a service call. Try the suggestions below to see if you can solve the problem before calling the service agent.

Problem	Possible Causes	Possible Solutions
The compressor does not start or make noise	The power is disconnected. The appliance is turned off. Blown fuse or circuit breaker Defective or broken wiring Overloaded outlet Defective relay	Check the power cord for proper connection. Turn the power ON. Replace the fuse or circuit breaker. Repair or replace wiring Replace the overloaded outlet Replace the relay
Warm storage temperatures	The appliance is in a location that is too warm. Low voltage, compressor cycling on overload	Relocate the appliance Check the voltage supply
The compressor runs continuously/ stored items are too cold	Defective control feeler tube. The control feeler tube was not installed properly	Call for service Push the control feeler tube well into the unit.
P1	Ambient probe failure	Call the service
P2	Evaporator probe failure	Call the service

If you continue to experience issues, contact Accucold's technical support department at **1-888-4-MEDLAB (1-888-463-3522)**.

NOTES

LIMITED WARRANTY

ONE-YEAR LIMITED WARRANTY

Within the 48 contiguous United States, for one year from the date of purchase, when this appliance is operated and maintained according to instructions attached to or furnished with the product, warrantor will pay for factory-specified parts and repair labor to correct defects in materials or workmanship. Service must be provided by a designated service company. Outside the 48 states, all parts are warranted for one year from manufacturing defects. Plastic parts, shelves and cabinets are warranted to be manufactured to commercially acceptable standards and are not covered from damage during handling or breakage.

5-YEAR COMPRESSOR WARRANTY

- 1. The compressor is covered for 5 years.
- 2. Replacement does not include labor.

ITEMS WARRANTOR WILL NOT PAY FOR:

- 1. Service calls to correct the installation of your appliance, to instruct you how to use your appliance, to replace or repair fuses or to correct wiring or plumbing.
- 2. Service calls to repair or replace appliance light bulbs or broken shelves. Consumable parts (such as filters) are excluded from warranty coverage.
- Damage resulting from accident, alteration, misuse, abuse, fire, flood, acts of God, improper installation, installation not in accordance with electrical or plumbing codes, or use of products not approved by warrantor.
- 4. Replacement parts or repair labor costs for units operated outside the United States.
- 5. Repairs to parts or systems resulting from unauthorized modifications made to the appliance.
- 6. The removal and reinstallation of your appliance if it is installed in an inaccessible location or is not installed in accordance with published installation instructions.

DISCLAIMER OF IMPLIED WARRANTIES - LIMITATION OF REMEDIES

CUSTOMER'S SOLE AND EXCLUSIVE REMEDY UNDER THIS LIMITED WARRANTY SHALL BE PRODUCT REPAIR AS PROVIDED HEREIN. IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO ONE YEAR. WARRANTOR SHALL NOT BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, OR LIMITATIONS ON THE DURATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, OR LIMITATIONS ON THE EXCLUSIONS OR LIMITATIONS MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE.

WARNING: This product can expose you to chemicals including Nickel (Metallic) which is known to the State of California to cause cancer.

For more information go to www.P65Warnings.ca.gov

Note: Nickel is a component in all stainless steel and some other metallic compositions.

Felix Storch, Inc. An ISO 9001:2015 registered company 770 Garrison Ave Bronx, New York 10474 www.accucold.com



For parts and accessory ordering, troubleshooting, and helpful hints, visit: www.summitappliance.com/support